

Central Venous Access & Dialysis



Central Venous Catheter Kit

CD - 1001

- Central Venous Catheter made of specially formulated and biocompatible Polyurethane material provides strength during insertion and also softens at body temperature to conform to the body tissues and reduces the risk and vascular trauma.
- Specially designed soft & beveled tip for smooth & easy insertion of catheter.
- Soft Flexible J-Tip Guide wire prevents the vessel perforation and also provides good torque to ensure film insertion.
- Sufficiently radio-Opaque material of catheter with clear, definite marking facilitates correct placement of catheter tip.



- Kink resistant Guidewire with soft & flexible J-tip offers better torque which helps in easy insertion & prevents vessel perforation.
- Wires are Double Distal.
- 2 piece design of guidewire advancer.
- Sterile / Disposable / Individually Tray Packed




Guidewire

Ordering Information

	Single Lumen Catheter	Catheter Size Ga / Fr	Lumen Ga	Dilator Fr	Guide Wire Nickel Titanium Double Distal	With Std. Needle	With Y Needle Std. Syringe	With Std. Needle Introducer Syringe
						P/N	P/N	P/N Catheter Size
Single Lumen Catheter								
1 Lumen 22 Ga 13 cms	22 Ga	22	13	4	0.018 x 45 cm	MCVC - 20 G 01 - 13	-	MCVC - 22 G 01 13 IS
1 Lumen 20 Ga 13 cms	20 Ga	20	13	4	0.018 x 45 cm	MCVC - 20 G 01 - 13	-	MCVC - 20 G 01 13 IS
1 Lumen 18 Ga 20 cms	18 Ga	18	20	5	0.018 x 45 cm	MCVC - 18 G 01 - 20	-	MCVC - 18 G 01 20 IS
1 Lumen 16 Ga 16 cms	16 Ga	16	16	6	0.035 x 45 cm	MCVC - 16 G 01 - 16	MCVC - 16 G 01 - 16 YN	MCVC - 16 G 01 16 IS
1 Lumen 16 Ga 20 cms	16 Ga	16	20	6	0.035 x 45 cm	MCVC - 16 G 01 - 20	MCVC - 16 G 01 - 20 YN	MCVC - 16 G 01 20 IS
1 Lumen 14 Ga 20 cms	14 Ga	14	20	8	0.035 x 45 cm	MCVC - 14 G 01 - 20	MCVC - 14 G 01 - 20 YN	MCVC - 14 G 01 20 IS
Double Lumen Catheter								
2 Lumen 4 Fr 6 cms	4 Fr	20, 22	6	5	0.018 x 45 cm	MCVC - 402 - 6	-	MCVC - 402 6 IS
2 Lumen 4 Fr 8 cms	4Fr	20, 22	8	5	0.018 x 45 cm	MCVC - 402 - 8	-	MCVC - 402 8 IS
2 Lumen 4 Fr 13 cms	4 Fr	20, 22	13	5	0.018 x 45 cm	MCVC - 402 - 13	-	MCVC - 402 13 IS
2 Lumen 5 Fr 8 cms	5 Fr	18, 20	8	6	0.018 x 45 cm	MCVC - 502 - 8	-	MCVC - 502 8IS
2 Lumen 5 Fr 13 cms	5 Fr	18, 20	13	6	0.018 x 45 cm	MCVC - 502 - 13	-	MCVC - 502 13 IS
2 Lumen 7 Fr 16 cms	7 Fr	14, 18	16	8	0.035 x 45 cm	MCVC - 702 - 16	MCVC - 702 - 16 YN	MCVC - 702 16 IS
2 Lumen 7 Fr 20 cms	7 Fr	14, 18	20	8	0.035 x 45 cm	MCVC - 702 - 20	MCVC - 702 - 20 YN	MCVC - 702 20 IS
Triple Lumen Catheter								
3 Lumen 4.5 Fr 6 cms	4.5 Fr	20, 22, 22	6	6	0.018 x 45 cm	MCVC - 4.503 - 6	-	MCVC - 4.503 - 6 IS
3 Lumen 4.5 Fr 10 cms	4.5 Fr	20, 22, 22	10	6	0.018 x 45 cm	MCVC - 4.503 - 10	-	MCVC - 4.503 - 10 IS
3 Lumen 5.5 Fr 6 cms	5.5 Fr	20, 22, 22	6	6	0.018 x 45 cm	MCVC - 5.503 - 6	-	MCVC - 5.503 - 6 IS
3 Lumen 5.5 Fr 8 cms	5.5 Fr	20, 22, 22	8	6	0.018 x 45 cm	MCVC - 5.503 - 8	-	MCVC - 5.503 - 8 IS
3 Lumen 5.5 Fr 13 cms	5.5 Fr	20, 22, 22	13	6	0.018 x 45 cm	MCVC - 5.503 - 13	-	MCVC - 5.503 - 13IS
3 Lumen 7 Fr 10 cms	7 Fr	16, 18, 18	10	8	0.035 x 45 cm	MCVC - 703 - 10	MCVC - 703 - 10 YN	MCVC - 703 - 10 IS
3 Lumen 7 Fr 13 cms	7 Fr	16, 18, 18	13	8	0.035 x 45 cm	MCVC - 703 - 13	MCVC - 703 - 10 YN	MCVC - 703 - 13 IS
3 Lumen 7 Fr 16 cms	7 Fr	16, 18, 18	16	8	0.035 x 45 cm	MCVC - 703 - 16	MCVC - 703 - 16 YN	MCVC - 703 - 16 IS
3 Lumen 7 Fr 20 cms	7 Fr	16, 18, 18	20	8	0.035 x 45 cm	MCVC - 703 - 20	MCVC - 703 - 20 YN	MCVC - 703 - 20 IS
Four Lumen Catheter								
4 Lumen 8.5 Fr 10 cms	8.5 Fr	14, 16, 18, 18	10	9	0.035 x 45 cm	MCVC - 8.504 - 10	MCVC - 8.504 - 10 YN	MCVC - 8.504 - 10 IS
4 Lumen 8.5 Fr 13 cms	8.5 Fr	14, 16, 18, 18	10	9	0.035 x 45 cm	MCVC - 8.504 - 13	MCVC - 8.504 - 13 YN	MCVC - 8.504 - 13 IS
4 Lumen 8.5 Fr 16 cms	8.5 Fr	14, 16, 18, 18	16	9	0.035 x 45 cm	MCVC - 8.504 - 16	MCVC - 8.504 - 16 YN	MCVC - 8.504 - 16 IS
4 Lumen 8.5 Fr 20 cms	8.5 Fr	14, 16, 18, 18	20	9	0.035 x 45 cm	MCVC - 8.504 - 20	MCVC - 8.504 - 20 YN	MCVC - 8.504 - 20 IS

Option of Needle-less valve is also available for Injection port

CVC kits are available with options of std. Needle, Needle with Back Flow valve & Introducer Syringe.



Complete Set of CVC Kit Consist

- Indwelling catheter.
- Y-Introducer needle with check valve.
- J-Tip guide wire.
- Vessel dilator.
- Luer lock syring.
- Scalpel.
- Catheter holder.
- Catheter holder clamp.
- Injection cap.
- Extension Line Clamp.
- Guiding Syringe / LOR Syringe (Optional).

