

## Hemodialysis Catheter Kit

CD - 1002

- Polyurethane material.
- Soft, geometrically designed conical tip to ensure easy insertion and prevent catheter related trauma.
- Hemodialysis catheter are single/double/multiple lumen catheter that provides temporary vascular access for hemodialysis until a permanent access is available or until another type of dialysis therapy is substituted.
- He multiple lumen catheter contains two large bore lumens that are connected to the dialysis machine to form a complete circuit for the removal and return of the patient's blood during treatment.
- Clear silicon lumen extensions for enhanced visibility and safety.  
Sterile / Disposable / Individually Tray Packed.



Part No	Dialysis Catheter	Dilator	Guidewire	Scalpel	Needle	Syringe
<b>Straight Extension</b>						
MDC - S107013 R	1 Lumen 7 F 13cms	8.5 F x 10 cm	0.035 x 70 cm	11 No. Thumb Blade	18 G 7 cms	5 ml Syringe
MDC - S108016 R	1 Lumen 8 F 16cms	9 F x 10 cm	0.035 x 70 cm	11 No. Thumb Blade	18 G 7 cms	5 ml Syringe
MDC - S206511 R	2 Lumen 6.5 F Hems	7 F x 10 cm	0.021 x 45 cm	11 No. Thumb Blade	18 G 7 cms	5 ml Syringe
MDC - S208511 R	2 Lumen 8.5 F Hems	9 F x 10 cm	0.021 x 45 cm	11 No. Thumb Blade	18 G 7 cms	5 ml Syringe
MDC - S211513 R	2 Lumen 11.5 F 13cms	11.5 F x 20 cm 10 F x 15 cms	0.035 x 75 cm	11 No. Thumb Blade	18 G 7 cms	5 ml Syringe
MDC - S212016 R	2 Lumen 12 F 16cms	11.5 F x 20 cm 10 F x 15 cms	0.035 x 75 cm	11 No. Thumb Blade	18 G 7 cms	5 ml Syringe
MDC - S312016 R	3 Lumen 12 F 16cms	11.5 F x 20 cm 10 F x 15 cms	0.035 x 75 cm	11 No. Thumb Blade	18 G 7 cms	5 ml Syringe
<b>Corved Extension</b>						
MDC - S206511 R	2 Lumen 6.5 F Hems	7 F x 10 cm	0.021 x 45 cm	11 No. Thumb Blade	18 G 7cms	5 ml Syringe
MDC - S208511 R	2 Lumen 8.5 F Hems	9 F x 10 cm	0.021 x 45 cm	11 No. Thumb Blade	18 G 7cms	5 ml Syringe
MDC - S211513 R	2 Lumen 11.5 F 13cms	11.5 F x 20 cm 10 F x 15 cms	0.035 x 70 cm	11 No. Thumb Blade	18 G 7cms	5 ml Syringe
MDC - S212016 R	2 Lumen 12 F 16cms	11.5 F x 20 cm 10 F x 15 cms	0.035 x 70 cm	11 No. Thumb Blade	18 G 7cms	5 ml Syringe
MDC - S312016 R	3 Lumen 12 F 16cms	11.5 F x 20 cm 10 F x 15 cms	0.035 x 70 cm	11 No. Thumb Blade	18 G 7cms	5 ml Syringe