




Blood Tubing Set

CD - 1003

- Medical grade PVC tubing for higher bio-compatibility.
- Advantage of built-in Heparin and Saline line.
- Extra corporeal blood circuit tubing set, used during dialysis
- Low Priming Volume.
- Easy of use.
- Customised blood tubing set as per clinic protocol and for different dialysis machines are available.

- Injection Site : Large finger guard.
- Drip Chamber : Available Various diameter sizes.
- Clamp : Clamps to secure blocking performance.
Sterile / Disposable / Individually Packed.



Peritoneal Dialysis Set

CD - 1005

- Super smooth kink resistance tubing ensures uniform flow rate.
- Clear, Transparent & Flexible built in Chamber with sharp and easy piercing Spike.
- Special designed for 'Y' shaped transfusion set for administrating dialysis solution.
- Provided with two upper control clamp to facilitate the attractive change of solution bottles.
- Top lower control clamps help for easy input and drainage of solutions.
Sterile / Disposable / Individual packed.

Peritoneal Transfusion Dialysis Kit

CD - 1004

- Suitable for performing peritoneal dialysis in patients of all age groups.
- Perforated open-end catheter manufactured from polypropylene.
- Stainless steel trocar is provided to facilitate smooth penetration.
- Scalpel blade is provided for the incision.
- The trocar and the tip of the catheter are perfectly matched to facilitate trauma free insertion.
- Latex flash bulb is provided on the junction unit for additional medication.
- Sizes : Adult, Child
Sterile / Disposable / Individually Tray Packed.



A V Fistula Needle

CD - 1006

[Fistula Needle with fixed wing,
Fistula Needles with rotating wing]

- Super smooth kink resistance tubing ensures uniform flow rate.
- To connect blood lines of the blood vessels through needles when dialysis is carried out through an internal fistula.
- Back eye needle to minimize interruption of blood flow.
- Siliconised ultra thin walled & sharp beveled needle to minimize trauma to the patient.
- Flexible butterfly wing for proper fixation.
- Color coded wing for easy identification of sizes as per standard.
- Also available in Twin pack.

Color Code	Blue	Green	Orange
Size G	15	16	17

