

## Rectal Catheter

U - 2024

- For introduction of enema solution into rectum to release/aspire rectal fluid.
  - Super smooth kink resistance tubing ensures uniform flowrate.
  - Atraumatic, soft rounded, closed tip with two lateral eyes for efficient drainage.
  - Frozen surface tubing for super smooth intubation.
  - Proximal end is fitted with universal funnel shaped connector for extension.
  - Color coded plain connector for easy identification of size.
  - Length : 40cm.
- Sterile / Disposable / Individually Packed.*


| Color Code | Red | Yellow | Violet | Dark Blue | White | Dark Green | Grey | Brown | Green | Light Green |
|------------|-----|--------|--------|-----------|-------|------------|------|-------|-------|-------------|
| Size FG    | 18  | 20     | 22     | 24        | 26    | 28         | 30   | 32    | 34    | 36          |




## Urethral Catheter

U-2025

- Used for short-term bladder catheterization through urethra.
  - Super smooth kink resistance tubing ensures uniform flowrate.
  - Atraumatic, soft rounded, closed tip with two lateral eyes for efficient drainage.
  - Frozen surface tubing for super smooth intubation.
  - Proximal end is fitted with universal funnel shaped connector for extension.
  - Color coded plain connector for easy identification of size.
  - X-Ray opaque line provided through out the length of catheter.
  - Length : 40cm.
- Sterile / Disposable / Individually Packed.*




| Color Code | Light Green | Blue | Black | White | Green | Orange | Red | Yellow | Violet | Light Blue |
|------------|-------------|------|-------|-------|-------|--------|-----|--------|--------|------------|
| Size FG    | 6           | 8    | 10    | 12    | 14    | 16     | 18  | 20     | 22     | 24         |



## Double J Stent

U - 2026

- D. J. Stent used for temporary internal drainage from ureteropelvic junction to the bladder
  - Made of superior grade poly-urethane material.
  - Black line on the stent indicates the direction of coil after withdrawal of guide wire .
  - Complete radio opaque stent for x-ray visualization. Implant tested with an externally low encrustation tendency.
  - Ureteral Stent Consist : Stent, Pusher, with or without Thread.
  - Ureteral Stent Kit Consist : Stent, Pusher, Guidewire, 2 Clamps with or without Thread.
- Sterile / Disposable / Individually Packed.*



|            |                              |
|------------|------------------------------|
| Size in FG | 4.0, 4.5, 4.8, 5.0, 5.5, 6.0 |
| Length CM  | 16, 20, 26, 28               |
| Type       | Open end & Close end         |