


Blood Transfusion Sets

T-1019 to 1023

- For transfusion of blood or blood components.
 - Super smooth kink resistance tubing ensures uniform flow rate.
 - Clear, Transparent & Flexible Drip Chamber with sharp and easy piercing Spike.
 - With 200 micron blood filter.
 - Smooth Roller Clamp Facilitates easy, safe control and adjustment of fluid rates.
- Sterile / Disposable / Individually Packed.




Non Vented, Double Chamber,
'Y' Connection & Luer Lock

Vented, Single Chamber,
Injection Bulb & Luer Lock



Non Vented, Double Chamber & Injections Bulb



OPTIONS AVAILABLE :

- Double Chamber or Single Chamber.
- With or Without Airvent.
- With or Without Needle.
- With or Without "Y" Connection.
- With or Without Luer Lock.
- Available with DEHP Free Material.
- Tube Length : 150cms, 180cms, 200cms.
- Size : Adult, Child.