

Oxygen Mask
A-3005



- Super smooth kink resistance tubing ensures uniform flowrate.
- Made from a soft clear plastic and are anatomically formed for comfort.
- Smooth and feathered edge of face mask for patient comfort and reducing irritation points.
- Adjustable elastic strap and integrated nose clip for proper positioning of mask.
- Proximal end of tube is connected with funnel shape connector for easy connection with oxygen source.
- Star lumen main tube to avoid accidental blockage.
- Tube Length : 210cm.
- Size : Adult, Paediatric & Neonatal.
- Available with DEHP free Material.
Sterile / Disposable / Individually Packed.



Guedel Airways

A-3006

- Suitable for maintaining an unobstructed oropharyngeal airway during general anaesthesia and in unconscious patients.
- Bite block to prevent biting of tongue and airway occlusion.
- Rounded atraumatic edges.
- Smooth airway path for easy cleaning.
- Color coded Bite Blocks for easy identification of sizes.
Sterile / Disposable / Individually Packed.

Color Code	Violet	Blue	Black	White	Green	Orange	Red	Yellow
Size	000	00	0	1	2	3	4	5



Suction Catheter With
Fingertip Control

A - 3008



- Super smooth kink resistance tubing ensures uniform flowrate.
- For removal of secretion from trachea and bronchial region.
- Atraumatic, soft rounded, Open Distal end with one lateral eyes for non-traumatic insertion.
- Frozen surface tubing for super smooth intubation.
- Proximal end is fitted with universal funnel shaped connector for extension.
- Color coded plain connector for easy identification of size.
- Available with DEHP free Material.
Sterile / Disposable / Individually Packed.

SUCTION CATHETER (PLAIN)

Provided with universal funnel shape connector for safe connection to standard suction equipment.

SUCTION CATHETER (Thumb Control with Stopper)

Suction Catheter with Thumb valve provided with stopper.

SUCTION CATHETER (Finger Tip)

Provided with Y shape vacuum control valve with finger tip control facility.

Color Code	Light Green	Blue	Black	White	Green	Orange	Red	Yellow	Violet	Light Blue
Size FG	6	8	10	12	14	16	18	20	22	24

Suction Catheter With
Thumb Control

A - 3007



Suction Catheter Plain

A - 3009

