

Multiflow Venturi Mask

A - 3017

- Super smooth kink resistance tubing ensures uniform flowrate.
- Smooth and feathered edge of face mask for patient comfort and reducing irritation points.
- Adjustable elastic strip and integrated nose clip for proper positioning of mask.
- Multiflow Venturi masks are devices that are constructed to supply oxygen or other gases to an individual.



- Masks fit snugly over the nose and mouth, and are equipped with 6 different oxygen concentration diluters which allows to select different oxygen concentration accordingly, and a tube that connects the oxygen mask to a storage tank where the oxygen is contained, to be connected to central oxygen supplying system in different hospitals.
- Proximal end of tube is connected with funnel shape connector for easy connection with oxygen source.
- Star lumen main tube to avoid accidental blockage.
- Tube Length : 210cm.
- Size : Adult, Paediatric & Neonatal.
- Available with DEHP free Material.
Sterile / Disposable / Individually Packed.

Venturi Color	Blue	Sky Blue	White	Yellow	Red	Orange	Green
Oxygen Concentration	24 %	28 %	31 %	35 %	40 %	50 %	60 %

Single Dial Venturi Mask

A - 3019

- Super smooth kink resistance tubing ensures uniform flowrate.
- Smooth and feathered edge of face mask for patient comfort and reducing irritation points.
- Adjustable elastic strip and integrated nose clip for proper positioning of mask.
- Single Dial Venturi masks are devices that are constructed to supply oxygen or other gases to an individual.
- There are two Oxygen diluters, one(White) for 24%, 26%, 28% and 30%, while the other (Blue) is for 35%, 40% or 50%. The appropriate percentage of Oxygen can be selected setting by the indicator of the diluter.
- Transparent (green or white) plastic masks also leave the face visible, allowing care providers to better ascertain patients' conditions.
- Proximal end of tube is connected with funnel shape connector for easy connection with oxygen source.
- Star lumen main tube to avoid accidental blockage.
- Tube Length : 210cm.
- Size : Adult & Paediatric.
- Available with DEHP free Material.
Sterile / Disposable / Individually Packed.



High Concentration Oxygen Mask

A - 3021

- Super smooth kink resistance tubing ensures uniform flowrate.
- Smooth and feathered edge of face mask for patient comfort and reducing irritation points.
- High concentration Oxygen mask with tubing is intended to supply oxygen or other gases to an individual for oxygen therapy, and the reservoir bag is able to avoid rebreathing.
- Proximal end of tube is connected with funnel shape connector for easy connection with oxygen source.
- Star lumen main tube to avoid accidental blockage.
- Tube Length : 210cm.
- Size : Adult & Paediatric.
- Available with DEHP free Material.
Sterile / Disposable / Individually Packed.



T Oxygenator

A - 3023

- T-Piece Kit Designed for use in recovery for the delivery of Oxygen via the patient 15mm ET Connection or Laryngeal mask.
- 7.6mm port for suction with non-popping cap.
- The kit comprises T-piece with port, 22mm – 15 cm reservoir flex tube and Oxygen stem with cap.
Sterile / Disposable / Individually Packed.

Options Available

- With or Without 2mtr. Star lumen Oxygen Tube.
- With or Without 40% or 60% Venturi Valve and 2 mtr. Tubing.

