

Suction Catheter Closed System (24 hr. / 72 hr.)

A - 3037

- Used to maintain ventilation and oxygen therapy while suctioning in mechanically ventilated patients.
- Minimises the risk of Ventilator Associated Pneumonia (VAP).
- Kit comes with 3 easily replaceable catheters, minimizing the need to replace entire suction system every 24-72 hrs.
- Catheter with lateral eyes, smooth tip and depth markings for maximum patient comfort.
- MDI port provides an effective, quick and convenient way of drug delivery.
- Irrigation port with one way valve allows lavage and irrigation of the suction catheter.
- On/off valve allows suction catheter tip to be isolated and creates an enclosed cleaning chamber.
- Transparent connectors and durable sleeve facilitates catheter manipulation and prevents infiltration of blood and sputum.
- 360° swivel elbow enables rotation for enhanced patient comfort and reduces torque.

Sterile / Disposable / Individually Packed.



Size	SIZES					
	6 Fr.	8 Fr.	10 Fr.	12 Fr.	14 Fr.	16 Fr.
O.D. (mm)	2.0	2.7	3.3	4.0	4.7	5.3l.
D. (mm)	1.4	1.8	2.4	2.8	3.2	3.8
Length (mm)	300	300	600	600	600	600
Color Code	Pea Green	Cambridge Blue	Black	White	Green	Orange



Epidural Needle

A - 3038

Quincke Point or Pencil Point
With or Without Introducer

- Highly suitable for Epidural Anaesthesia.
- Markings on the needle at every 10 mm distance from tip making it easy to observe and also control the exact depth of insertion. Tuohy type needle tip with fixed wings.
- Clear polycarbonate hub enables visualization of Cerebro Spinal Fluid (CSF) or blood easily.
- Polished inner bevel minimizes the risk of catheter shearing.
- Standard Length : 80mm.
- Sterile / Disposable / Individually Packed.

Color Code	Blue	Pink	Bige
Size G	16	18	19

Epidural Catheter

A - 3039

- Epidural Catheter is specially designed for short term & long term continuous anaesthesia & pain therapy.
 - Catheter is made from unique plasticizer – free, latex – free formulation of polymer with excellent bio compatibility.
 - Flexible & atraumatic soft tip with three lateral eyes reduces the risk of complications.
 - Catheter is printed to determine the depth of insertion & accurate placement of catheter in the Epidural space.
 - Crystal clear catheter with radio opaque line provides clear visualization of blood / CSF.
 - Luer lock twist connector is supplied for easy & safe connection.
 - Tube Length : 80 cms.
 - Sizes : Passes through G 16, 18 needle.
- Sterile / Disposable / Individually Packed.

