



Twin Bore Nasal Oxygen Set

A-3001 - 3002


- Super smooth kink resistance tubing ensures uniform flowrate.
 - Designed for easy administration of oxygen and comfort of patient.
 - Soft twin prong nasal tips to ensure equal volume of oxygen to both air passages.
 - Star lumen main tube to avoid accidental blockage.
 - Smoothly finished and adjustable nasal tips for maximum patient comfort.
 - Soft funnel shaped connector facilitates easy connection to oxygen source.
 - Tube Length : 210 cm.
 - Size : Adult, Paediatric & Neonatal.
 - Available with DEHP free Material.
- Sterile / Disposable / Individually Packed.



Nasal Oxygen Catheter

A-3003

- Super smooth kink resistance tubing ensures uniform flowrate.
 - Suitable for direct administration of oxygen via nasopharyngeal route Atraumatic, soft rounded, Open Distal end with two lateral eyes for non-traumatic insertion and dispersion of oxygen and prevention of oxygen burn.
 - Frozen surface tubing for super smooth intubation.
 - Proximal end is fitted with universal funnel shaped connector for connection to oxygen source.
 - Color coded plain connector for easy identification of size.
 - Available with DEHP free Material.
- Sterile / Disposable / Individually Packed.



Color Code	Light Green	Blue	Black	White	Green	Orange	Red	Yellow	Violet	Light Blue
Size FG	6	8	10	12	14	16	18	20	22	24