




Ryle's Tube

G - 5001

- Suitable for nasogastric insertion for nutritional purposes or aspiration of intestinal secretion.
 - Super smooth kink resistance tubing ensures uniform flowrate.
 - Atraumatic, soft rounded, closed tip with four lateral eyes for efficient aspiration & administration.
 - The tube is marked at 50,60 & 70cm from the tip for accurate placement.
 - Proximal end is fitted with universal funnel shaped connector for extension.
 - Color coded connectors for easy identification of size as per standards.
 - X-Ray opaque line provided through out the length of Tube.
 - Available with DEHP Free Material.
 - Available with universal luer Connector.
 - Tube Length : 105cm.
- Sterile / Disposable / Individual Packed.



Color Code	Light Green	Blue	Black	White	Green	Orange	Red	Yellow	Violet	Light Blue
Size FG	6	8	10	12	14	16	18	20	22	24

Stomach Tube

G - 5002

- Super smooth kink resistance tubing ensures uniform flowrate.
 - Atraumatic, soft rounded, closed tip with four lateral eyes for efficient aspiration & administration.
 - Tube is marked at 45, 55, 65 and 75 cm from the tip for accurate placement.
 - Proximal end is fitted with universal funnel shaped connector for extension.
 - Color coded connectors for easy identification of size as per standards.
 - X-Ray opaque line provided through out the length of Tube.
 - Available with DEHP Free Material.
 - Available with universal luer Connector.
 - Open Distal end or closed distal end.
 - Tube Length : 105cm.
- Sterile / Disposable / Individual Packed.



Color Code	Light Green	Blue	Black	White	Green	Orange	Red	Yellow	Violet	Light Blue
Size FG	6	8	10	12	14	16	18	20	22	24