



Karman Cannula

M - 6001 - 6002

- Specially designed for aseptic medical termination of pregnancy.
 - Distal end is coned with two large lateral eyes to facilitates the curette.
 - Suitable for use with MTP syringe or suction apparatus.
 - With or Without universal adopter for connecting to suction apparatus.
- Sterile / Disposable / Individual Packed.




Color Code	Grey	Orange	Red	Green	Blue	Light Blue
Size MM	5	6	7	8	10	12

Umbilical Cord Clamp

M - 6003

- Suitable for clamping the umbilical cord of new born baby immediately after the birth.
- Provided with double purpose security lock click to indicate the correct locking and protect against accidental reopening.
- Finger grip ensures safe and convenient handling, particularly when gloves are wet.
- Provided with grooves all along the length to prevent the slipping of the umbilical cord and to retain it in the same position.
- Supplied in Blue and White colors, other colors are available upon demand.

Sterile / Disposable / Individually Packed.



Corrugated Drainage Sheet

M - 6005

- Extra soft super smooth PVC corrugated sheet for wound drainage.
- X-ray opaque line provided throughout the length of sheet.
- Multi channel drainage system of efficient drainage during operation.

Sterile / Disposable / Individually Packed.

