

**Examination Polyethylene Gloves**  
M - 6006 - 6007



- Made of 100% Virgin high density polyethylene.
- High tearing strength.
- Safe & Hygienic.

**Latex Surgical Gloves**  
M - 6008, M - 6009

- Surgical Gloves are designed to be perfectly anatomical for comfortable wrinkle-free fitting.
  - The extremely soft rubber reduces finger fatigue and increases efficiency.
  - The special rubber film used in every pair of Gloves adds sensitivity while giving complete protection to the hand.
  - Apart from ergonomic features, our stringent quality control measures make it the safest protective hand wear for surgeons.
- Sterile or Non-Sterile**

**Special Features**

- Lower modulus for better comfort
- Low powder content that minimizes skin irritation
- Thicker gloves for enhanced safety

Size	6	6 ½	7	7 ½	8
------	---	-----	---	-----	---



**Examination Rubber Gloves**  
M - 6010

- 100% natural rubber latex.
- Convenient and comfortable to wear.
- Non finger fatigue.
- Keeps the hands cool and comfortable.
- Examination gloves are available in 3 materials Latex, Nitrile (Latex Free) & Vinyl (Latex Free) type. The examination gloves are tested and proved safe against dermal sensitization and skin irritation.
- Pre-powdered or Powder Free

Size	Small, Medium, Large, Extra Large
------	-----------------------------------



**Disposable Face Mask**  
M - 6013 - 6014

- Manufactured from soft non-woven cloth.
  - Double layer mask offering excellent bacterial filtration efficiency.
  - Super smooth soft ultrasonically welded mask is most comfortable during use.
  - Tie On or Ear loop type.
- Sterile or Non-Sterile.**

**Infant Mucus Extractor**  
M - 6015

- Specially designed for aspiration of secretion from Oropharynx in newly born babies to ensure free aspiration.
  - Clear transparent container permits immediate visual examination of the aspirate.
  - Also suitable for obtaining mucus specimen for micro biological examination.
  - Spare closure cap is provided to seal the container for safe transportation of specimen to the laboratory or aseptic disposal of container.
  - Low friction surface catheter is provided with open end silk smooth round tip for trauma free insertion.
  - Capacity 20 to 25ml.
- Sterile / Disposable / Individually Packed**

Catheter	Size	Length
Standard	10 FG	38cm
Optional	12 FG, 14 FG	50cm



**Disposable Cap**  
M - 6011 - 6012



- Manufactured from soft non-woven cloth.
  - Specially designed to protect the patient against dropping of hair / dandruff of the doctors and other attending staff.
  - Provided with stitched interlocking adaptable elastic for better grip on the forehead and comfortable fitting.
- Sterile or Non-Sterile.**