

### Blood Donor Sets

T-1025



- Super smooth kink resistance tubing ensures uniform flow rate.
- Smooth Roller Clamp Facilitates easy, safe control and adjustment of fluid rates.
- Needle size G: 16. extra sharp needle for smooth and effortless vein puncture.  
Sterile / Disposable / Individually Packed.



### Insuline Syringe

T-1043

- Ideal sharpness for less pain and more comfort
- Advanced Uni-Body design.
- Low dead space design ensures dosage accuracy and minimum wastage of drugs.
- Clear transparent barrel and bold marking facilitates proper drug delivery.
- Needles have ultra thin wall permitting rapid and efficient delivery.  
Sterile / Disposable / Individually Packed

Nominal Capacity	Strength of Insulin	Needle Gauge	Needle Length	Color Code of Top & Bottom Cap
1 ml	U - 40	29 G	13 mm	Red
1 ml	U - 40	30 G	8 mm	Red
1 ml	U - 40	31 G	8 mm	Red
1 ml	U - 100	29 G	13 mm	Orange
1 ml	U - 100	30 G	8 mm	Orange
1 ml	U - 100	31 G	8 mm	Orange





### Disposable Syringe

T-1031 to 1042

[ With or Without Needle ] [ Luer Slip / Luer Lock ]

- Sufficiently transparent barrel allows easy measurement of the volume contained in the disposable syringe friction free.
- Thermoplastic Elastomer (TPE) is used to manufacture the Gasket for smooth movement of plunger.
- In superior quality by high accurate moulds.
- Needles are manufactured by using imported Cannula having precise tip sharpness and good bonding strength with the hub.
- Needles have ultra thin wall permitting rapid and efficient drug delivery.  
Sterile / Disposable / Individually Packed.

Size ML.						
1	2	3	5	10	20	50



Luer Lock Syringe

### IV Cannula T-1045

- The flash back chamber allows quick visualization of venous return.
- Wings specially designed for easy gripping and safe clamping.
- Needle cover to prevent accidental damage to needle & catheter.
- Color coded body/port cap as per ISO Standard.
- 6% luer taper as per International standard.

- Japanese stainless steel (AISI 304) needle with diameter and length adhering to all ISO standards.
- Triple point Needle for painless insertion.
- Thin walled, Radio-opaque stripped / transparent Teflon catheter, as per USP standard Class VI.
- Customized Automated Tipping Technology for lower penetration forces.
- Optimum trim distance for pain free penetration.  
Sterile / Disposable / Individually Packed.



With Port & With Wings



Without Port & Without Wings

Size	Color	Catheter	Size	Length	Flow Rate
14 G	Orange	1.7	2.1	45	270
16 G	Grey	1.3	1.7	45	180
17 G	White	1.1	1.5	45	125
18 G	Green	0.9	1.3	45	80
20 G	Pink	0.8	1.1	32	54
22 G	Blue	0.6	0.9	25	33
24 G	Yellow	0.5	0.7	19	20
26 G	Violet	0.4	0.6	19	15



Without Port & With Wings