




Close Wound Suction Unit

S - 4001

[MEDI-VAC SET]

- Suitable to offer surgeons and doctors an effective device for closed wound drainage under negative pressure post operatively with the option to use one catheter or two catheters simultaneously.
- Super smooth kink resistance tubing enables uniform flowrate.
- Radon drains catheters are provided with radio opaque line and satin smooth eyes.
- Connecting tube is kink resistant and is provided with additional strength to withstand the suction.
- Easy to use by one person depress the chamber to active the suction of bellow unit.
- Available in different sizes with matching size curved needle to meet moderate to heavy drainage needs.
- X-ray opaque line provided throughout the length of Tube.
- Connecting tube with clamp & "Y" connector, curved needle with matching catheter & spare perforated catheter.
- Available with DEHP Free Material.
Sterile / Disposable / Individually Packed.



Color Code	Black	White	Green	Orange	Red
Size FG	10	12	14	16	18

Mini Close Wound Suction Unit

S - 4002

BABY-VAC SET

- Baby Vac Set has been specially designed for Wound Drainage under negative pressure for minor surgery and plastic surgery.
- Super smooth kink resistance tubing ensures uniform flowrate.
- Easy to use Catheter to connect bellow container.
- Catheters are available in two sizes with matching size curved needle .
- Frozen surface Catheter and smooth round eyes facilitate efficient drainage.
- X-ray opaque line provided throughout the length of Tube.
- Curved Needle with matching catheter.
- Available with DEHP Free Material.
Sterile / Disposable / Individually Packed.



Color Code	Green	Red
Size FG	06	08