

Chest Drainage Catheter

S - 4003 - 4004

Thoracic Catheter [Straight or Curved]
[With or Without Trocar]

- Specially designed for rapid non operative pleural and chest drainage.
- Super smooth kink resistance tubing ensures uniform flowrate.
- Catheter without Trocar : Atraumatic, soft rounded, Open Distal end with six lateral eyes for non-traumatic insertion.
- Catheter With Trocar : Atraumatic, soft rounded, Open Distal end with two lateral eyes for non-traumatic insertion.
- Cross side eyes to prevent tissue aspiration.
- Large smooth drainage eyes for efficient drainage.
- Proximal end is fitted with tapered connector for easy connection to drainage bottle.
- X-ray opaque line provided throughout the length of catheter.
- Marking on catheter is provided at from the distal end to help ascertain.
- Matching size loose connector which facilitates easy connection to drainage system.
- Available with DEHP Free Material.

Size FG

08, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40



Chest Drainage Bottle

S - 4005

- Chest Drainage Bottle is a single use disposable unit which is designed to forward the collected Drainage into a collection chamber by the liquid static suction effects.
- Super smooth kink resistance tubing ensures uniform flowrate.
- Makes it easy to observe and measure the collected drainage content during surgery.
- Facilitates safe Retrograde Technique.
- Product is suitable for pleural drainage in conjunction with chest drainage catheter in Cardiac and Thoracic procedures.
- Easy to read graduation helps to determine the drain volume precisely.
- Product Is provided with hanging at a higher position as compares to the common collection chamber allowing the collection of drainage without giving rise to any inconvenience.
- Suction port is provided to connect with suction unit.
- Easy to read graduation helps to determine the drain volume precisely.
- Available with DEHP Free Material.

Sterile / Disposable / Individually Packed.

Size ML.

Single Chamber : 500, 1200, 1600, 2000
Double Chamber : 2500



Chest Drainage Kit

S - 4006

- Kit Includes all the essentials to perform chest drainage and permit the safe mobilisation of patients when under water seals cannot be established or are not practical.
 - Super smooth kink resistance tubing ensures uniform flowrate.
 - Ambulatory Chest Drainage System with integral tubing.
 - Available with DEHP Free Material.
- Sterile / Disposable / Individually Packed.

Abdominal Drain Kit

S - 4008 [Medi - ADK]

- Specially designed for post operative abdominal drainage.
 - Super smooth kink resistance tubing ensures uniform flowrate.
 - Soft abdominal drainage catheter with collection bag of 2000 ml. capacity.
 - Catheter is Atraumatic, soft rounded, Open Distal end with six lateral eyes for non-traumatic insertion.
 - X-ray opaque line provided throughout the length of catheter.
 - Specially designed handle holds tube up right and facilitates carrying.
 - Transparent sheeting allows visual inspection.
 - Available with DEHP Free Material.
- Sterile / Disposable / Individually Packed.

Size FG

16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36

